

# Disposable Endoscopic Hemoclip



## Disposable Endoscopic Hemoclip



### Intended use

- Compatible with endoscope, the device is indicated for clip placement within digestive tract for the purpose of mechanical pressure treatment of bleeding of small arteries and pulsation.

### Features

- Rotatable clips design, offering a precise placement.
- Repeated opening and closing clips.
- Innovative design of handle lock system and atraumatic tip, avoid damage to scope after releasing clips.
- Sensitive release system, easy to release the clips.



### Specifications

Model	O.D. (mm)	WorkingLength (mm)	Clip opening (mm)	Endoscope Channel (mm)
JHY-CD11-25-160	2.5	1600	11	≥2.8
JHY-CD11-25-180	2.5	1800	11	≥2.8
JHY-CD11-25-230	2.5	2300	11	≥2.8
JHY-CD13-25-160	2.5	1600	13	≥2.8
JHY-CD13-25-180	2.5	1800	13	≥2.8
JHY-CD13-25-230	2.5	2300	13	≥2.8
JHY-CD16-25-160	2.5	1600	16	≥2.8
JHY-CD16-25-180	2.5	1800	16	≥2.8
JHY-CD16-25-230	2.5	2300	16	≥2.8