

## Bite Block



## Biopsy Valve



## Polyp Trap



### Specifications

#### Bite Block

Model	Characteristics
JHY-KD-A1	Disposable, without strap, suitable for adult
JHY-KD-B1	Disposable, with strap, suitable for adult
JHY-KD-C1	Disposable, with strap, suitable for pediatric

#### Biopsy Valve

Model	Characteristics
JHY-BV-A1	Disposable, for Olympus & Fujinon Endoscope
JHY-BV-B1	Disposable, for Pentax Endoscope

#### Polyp Trap

Model	Characteristics
JHY-PT-05	Disposable, 5 chambers

Changzhou JiuHong Medical Instrument Co.,Ltd.