

Disposable Cytology Brush

Features

- Innovational brush design, with no risk of drop-off
- Straight-shaped brush: easy to enter the depths of respiratory and digestive tract with
- U-shaped brush: obtain tissues from all directions with a large quantity



Specifications

Intended use

- Tissue sampling in digestive and respiratory tracts

Model	Catheter O.D. (mm)	Endoscope Channel (mm)	Working Length (mm)	Brush O.D. (mm)	Brush length (mm)	Characteristics
JHY-BC-18-105-A3	1.8	≥2.0	1050	3	10	A: Straight Shaped Brush
JHY-BC-18-120-A3	1.8	≥2.0	1200	3	10	
JHY-BC-23-160-A3	2.3	≥2.8	1600	6	8	
JHY-BC-23-230-A3	2.3	≥2.8	2300	6	8	
JHY-BC-18-105-B	1.8	≥2.0	1050	3	10	B: "U" Shaped Brush
JHY-BC-18-120-B	1.8	≥2.0	1200	3	10	
JHY-BC-23-160-B	2.3	≥2.8	1600	6	8	
JHY-BC-23-230-B	2.3	≥2.8	2300	6	8	

A **Authorised Representative**
 IBS MEDICAL EQUIPMENTS SDN BHD
 Unit 07-05, Level 07, Tower A, Vertical Business Suite,
 Avenue 3, Bangsar South, No. 8, Jalan Kerinchi,
 59200, Kuala Lumpur, Malaysia.

T Tel : +6 03 2242 4024

W Whatsapp : +6 012 657 2400

W Website : www.ibsmedic.com

E Email : sales@ibsmedic.com