




Video Bronchoscope B41

MDA Reg. No. : GB 5960121-82276



Outer Diameter	Instrument Channel
3.2 mm	1.2 mm
4.0 mm	1.8 mm
5.2 mm	2.6 mm



Authorised Representative :

IBS Medical Equipments Sdn Bhd
Unit 07 - 05, Level 07, Tower A, Vertical Business Suite
Avenue 3, Bangsar South No.8, Jalan Kerinchi 59200
Kuala Lumpur, Malaysia.

T : +603 2242 3024  +6012 657 2400

www.ibsmedic.com
sales@ibsmedic.com

- Our scope is used for endoscopy and endoscopic surgery within the airways & tracheobronchial tree.
- Portable design for daily use
- With new features: chip-on-tip
- Totally no optical - fibers
- Extremely flexible & durable
- Long service life & low cost maintenance
- Anti-foggy: no need to pre-heat
- Image sensor and two LED lights on the distal end
- White Balance Control
- 5 Level Brightness Control
- Video / Image Recording
- Connection to HDMI Monitor



Purchase	Instalment (12 Months)	Rental (Monthly)
<ul style="list-style-type: none"> • 1 Year Warranty • 3 Weeks Delivery Upon PO 	<ul style="list-style-type: none"> • 1 Year Warranty • Free loaner unit within 48 working hours, returned once completed 	<ul style="list-style-type: none"> • Complete Set Rental • Replacements within 48 working hours



* 1 Year Warranty start from Shipping Date from Manufacturer. Estimated 2 weeks for delivery.

* Warranty covers manufacturing defects only.

<ul style="list-style-type: none"> • Experience and Professional Team • Ready Spare Parts 	<ul style="list-style-type: none"> • Fast Repair Time • Lowest Cost of Maintenance
---	--



Specification

Endoscope	Monitor \ Image Processor	
Brand : MDH Model : B41 Working Length : 600 mm Angulation : U180° / D130° Recording : Yes Resolution : 400*400 Outer Diameter : 3.2 / 4.0 / 5.2 mm Instrument Channel : 1.2 / 1.8 / 2.6 mm	 4.3" Monitor	 10.1" Monitor
	Storage : 8GB Recording : Yes Resolution : 800 x 400 (RGB) Power : 3.7V lithium battery, rechargeable Operating Time : 3 hours	Storage : 32GB Recording : Yes Resolution : 1280 x 800 (RGB) Power : 100-240V 50/60Hz Operating Time : >4 hours

MOF : K22184098821841368

MDA : KP11442186317

GDPMD MASFMD : 101035



Authorised Representative :

IBS Medical Equipments Sdn Bhd
 Unit 07 - 05, Level 07, Tower A, Vertical Business Suite
 Avenue 3, Bangsar South No.8, Jalan Kerinchi 59200
 Kuala Lumpur, Malaysia.

T : +603 2242 3024  +6012 657 2400

www.ibsmedic.com
sales@ibsmedic.com