

# Single-Use SUV-2A-B Urology Videoscope



**Improve image contrast**  
Gain off    Gain on

**White Balance**

**Edge structure and capillary emphasis, better identify upper tract tumors due to improved visibility of vessel structures.**  
S1    S2

**Photo**

**Video**

**Freeze image**

**Electronic Magnification**

**Three levels of adjustable brightness**

**USB storage**

**Scivita brings a brand-new and higher-resolution digital image without distortion, even on a 32-inch monitor.**

The processor compatible with various monitors. 32-inch monitor is highly recommended.

## SPECIFICATION

## SUV-2A-B



Mould	SUV-2A-B
Shaft Diameter OD	8.4 Fr(Φ2.8mm)
Outer Diameter of distal End	7.8 Fr
Weight	200g
Working Length	680mm
Total Length	990mm
Working Channel OD	3.6 Fr(Φ1.2mm)
Field of View	120°
Depth of Field	0.5-120mm
Angulation Range	285° Up/Down

## SPECIFICATION

### PROCESSOR



Mould	HDVS-S100A
Size	300mm x 225mm X 57mm
Weight	2.5 Kg
Power Supply Requirement	100-240V
Output	DVI*1、SDI*2
Interface	Foot Switch and USB


### ACCESSORIES(OPTIONAL):

Irrigation Valve	80460	39mm
	T-001	41mm
	Y-001	84mm
Foot Pedal	KF2 1S-MED GP25	



# Single-Use **SUV-2A-B** Urology Videoscope

MDA : GB10740222-103413



## Be ready at any time

- Plug and play. Off the shelf ready. Consistency every time.
- You decide which hand drives.
- Avoid cross-infection.
- Save on re-sterilisation and maintenance costs.
- Eradicate unnecessary repair bills.
- Save time post procedure, no need for scope integrity testing.



MOF : K22184098821841368



MDA : KP11442186317



GDPMD MASFMD : 101035

