




Video Flexible Cystoscope C41

MDA Reg. No. : GB 5723521-82280



Outer Diameter	Instrument Channel
5.2 mm	2.6 mm




Authorised Representative :

IBS Medical Equipments Sdn Bhd
No. 7-G, 3rd Mile Square Business Centre,
No. 151, Old Klang Road, 58100 Wilayah Persekutuan
Kuala Lumpur.

T : +603 7971 4221  +6012 346 3400

www.ibsmedic.com
sales@ibsmedic.com

- Our scope is designed for endoscopic access to and examination of the lower urinary tract including the bladder
- Portable design for daily use
- With new features: chip-on-tip
- Totally no optical - fibers
- Extremely flexible & durable
- Long service life
- Low cost maintenance
- Anti-foggy: no need to pre-heat
- Image sensor and two LED lights on the distal end
- White Balance Control
- 5 Level Brightness Control
- Video / Image Recording
- Connection to HDMI Monitor



Purchase	Instalment (12 Months)	Rental (Monthly)
<ul style="list-style-type: none"> • 1 Year Warranty • 3 Weeks Delivery Upon PO 	<ul style="list-style-type: none"> • 1 Year Warranty • Free loaner unit within 48 working hours, returned once completed 	<ul style="list-style-type: none"> • Complete Set Rental • Replacements within 48 working hours



* 1 Year Warranty start from Shipping Date from Manufacturer. Estimated 2 weeks for delivery.

* Warranty covers manufacturing defects only.

<ul style="list-style-type: none"> • Experience and Professional Team • Ready Spare Parts 	<ul style="list-style-type: none"> • Fast Repair Time • Lowest Cost of Maintenance
---	--



Specification

Endoscope	Monitor \ Image Processor	
Brand : MDH Model : C41 Working Length : 380 mm Angulation : U210° / D130° Recording : Yes Resolution : 400*400 Outer Diameter : 5.2 mm Instrument Channel : 2.6 mm	 4.3" Monitor	 10.1" Monitor
	Storage : 8GB Recording : Yes Resolution : 800 x 400 (RGB) Power : 3.7V lithium battery, rechargeable Operating Time : 3 hours	Storage : 32GB Recording : Yes Resolution : 1280 x 800 (RGB) Power : 100-240V 50/60Hz Operating Time : >4 hours

MOF : K22184098821841368

MDA : KP11442186317

GDPMD MASFMD : 101035



Authorised Representative :

IBS Medical Equipments Sdn Bhd
No. 7-G, 3rd Mile Square Business Centre,
No. 151, Old Klang Road, 58100 Wilayah Persekutuan
Kuala Lumpur.

T : +603 7971 4221  +6012 346 3400

www.ibsmedic.com
sales@ibsmedic.com