




Video Portable ENT Endoscope E41

MDA Reg. No. : GB 5739521-82276



Outer Diameter	Instrument Channel
2.8 mm	None
3.8 mm	None




Authorised Representative :

IBS Medical Equipments Sdn Bhd
Unit 07 - 05, Level 07, Tower A, Vertical Business Suite
Avenue 3, Bangsar South No.8, Jalan Kerinchi 59200
Kuala Lumpur, Malaysia.

T : +603 2242 3024  +6012 657 2400

www.ibsmedic.com
sales@ibsmedic.com

- Our Video Portable ENT Endoscope has been designed to be used with a control section, insertion tube and LCD Screen (including 4.3" or 10.1" two sizes for choosing), accessories and other ancillary equipment for diagnosis and treatment of diseases in ear, nose and throat.
- White balance control.
- 5 Level of brightness control.
- Recording of (images and videos in USB).
- Connection with HDMI Monitor.
- Portable designed for daily use.
- With new features: chip-on-tip.
- Totally no optical – fibres.
- Extremely flexible & durable.
- Long service life with low cost maintenance.
- Anti-foggy (No need to pre-heat).
- Image sensor and two LED lights on the distal end.



Purchase	Instalment (12 Months)	Rental (Monthly)
<ul style="list-style-type: none"> • 1 Year Warranty • 3 Weeks Delivery Upon PO 	<ul style="list-style-type: none"> • 1 Year Warranty • Free loaner unit within 48 working hours, returned once completed 	<ul style="list-style-type: none"> • Complete Set Rental • Replacements within 48 working hours



* 1 Year Warranty start from Shipping Date from Manufacturer. Estimated 2 weeks for delivery.

* Warranty covers manufacturing defects only.

<ul style="list-style-type: none"> • Experience and Professional Team • Ready Spare Parts 	<ul style="list-style-type: none"> • Fast Repair Time • Lowest Cost of Maintenance
---	--



Specification

Endoscope	Monitor \ Image Processor	
Brand : MDH Model : E41 Working Length : 330 mm Angulation : U180°/D130° Recording : Yes Resolution : 400*400 Outer Diameter : 2.8 & 3.8 mm Instrument Channel : Without Channel	 4.3" Monitor	 10.1" Monitor
	Storage : 8GB Recording : Yes Resolution : 800 x 400 (RGB) Power : 3.7V lithium battery, rechargeable Operating Time : 3 hours	Storage : 32GB Recording : Yes Resolution : 1280 x 800 (RGB) Power : 100-240V 50/60Hz Operating Time : >4 hours

MOF : K22184098821841368

MDA : KP11442186317

GDPMD MASFMD : 101035



Authorised Representative :

IBS Medical Equipments Sdn Bhd
 Unit 07 - 05, Level 07, Tower A, Vertical Business Suite
 Avenue 3, Bangsar South No.8, Jalan Kerinchi 59200
 Kuala Lumpur, Malaysia.

T : +603 2242 3024  +6012 657 2400

www.ibsmedic.com
sales@ibsmedic.com