




Video Flexible Ureteroscope U41

MDA Reg. No. : GB 5200122-83576



Outer Diameter	Instrument Channel
3.2 mm	1.2 mm

Authorised Representative :

IBS Medical Equipments Sdn Bhd
No. 7-G, 3rd Mile Square Business Centre,
No. 151, Old Klang Road, 58100 Wilayah Persekutuan
Kuala Lumpur.

T : +603 7971 4221  +6012 346 3400



www.ibsmedic.com
sales@ibsmedic.com

- Has been designed to be used with a control section insertion tube and LCD Screen, accessories and other ancillary equipment for visualization purposes during diagnostic and therapeutic procedures.
- Portable design for daily use
- With new features: chip-on-tip
- Totally no optical - fibers
- Extremely flexible & durable
- Long service life & low cost maintenance
- Anti-foggy: no need to pre-heat
- Image sensor and two LED lights on the distal end
- White Balance Control
- 5 Level Brightness Control
- Video / Image Recording
- Connection to HDMI Monitor



Purchase

- 1 Year Warranty
- 3 Weeks Delivery Upon PO

* 1 Year Warranty start from Shipping Date from Manufacturer. Estimated 2 weeks for delivery.
 * Warranty covers manufacturing defects only.

- Experience and Professional Team
- Ready Spare Parts

- Fast Repair Time
- Lowest Cost of Maintenance



Specification

Endoscope

Brand : MDH
 Model : U41
 Working Length : 685 mm
 Angulation : U275° /D275°
 Recording : Yes
 Resolution : 400*400
 Outer Diameter : 3.2 mm
 Instrument Channel : 1.2 mm

Monitor \ Image Processor



4.3" Monitor

Storage : 8GB
 Recording : Yes
 Resolution : 800 x 400 (RGB)
 Power : 3.7V lithium battery, rechargeable
 Operating Time : 3 hours



10.1" Monitor

Storage : 32GB
 Recording : Yes
 Resolution : 1280 x 800 (RGB)
 Power : 100-240V 50/60Hz
 Operating Time : >4 hours

MOF : K22184098821841368

MDA : KP11442186317

GDPMD MASFMD : 101035



Authorised Representative :

IBS Medical Equipments Sdn Bhd
 No. 7-G, 3rd Mile Square Business Centre,
 No. 151, Old Klang Road, 58100 Wilayah Persekutuan
 Kuala Lumpur.

T : +603 7971 4221 WhatsApp +6012 346 3400

www.ibsmedic.com
sales@ibsmedic.com