

Irrigation Pump

MDA Reg No : GB6159723-142636

Used for liquid irrigation in Bladder or Ureteric Urology surgery, capable to be used in Arthroscopic, Laparoscopic, Choledochoscopic and Hysteroscopic surgery.

Accessories

- JG100-01 (Irrigation Tube)
- JG100-05 (Connector)



Technical Parameters

- Supply voltage : ~ 220V 50 Hz
- Power : 150VA
- Pressure Range : 50 - 400mmHg (6.6 ~ 523.3kPa)
- Flowrate Range : 0.1 - 1.0L/min
- Noise: ≤ 70 dB (A)
- Brand : Jingrui

Authorised Representative :

IBS Medical Equipments Sdn Bhd
No. 7-G, 3rd Mile Square Business Centre, No. 151,
Old Klang Road, 58100 Wilayah Persekutuan
Kuala Lumpur.

T : +603 7971 4221

+6012 346 3400



Functional Characteristic

- Usage of roller pump for a larger flow and irrigation, tubings can be autoclaved to prevent secondary contamination.
- Pressure is adjustable on request and monitored in real time.



www.ibsmedic.com

sales@ibsmedic.com