


# Disposable Biopsy Forceps With Alligator Teeth



## Disposable Biopsy Forceps With Alligator Teeth



### Intended use

- Tissue sampling in digestive and respiratory tracts

### Features

- Alligator teeth design with a larger sampling capacity, no slipping during sampling.
- Easier to sample tissues, effectively avoid tearing.
- Smooth jaws, effectively decrease the friction to working channel.



### Specifications

Model	O.D. (mm)	Length (mm)	Endoscope Channel (mm)	Characteristics
JHY-FB-18-105-O-O-A1	1.8	1050	≥ 2.0	without spike, uncoated, with alligator teeth
JHY-FB-18-105-O-P-A1	1.8	1050	≥ 2.0	without spike, coated, with alligator teeth
JHY-FB-18-160-O-O-A1	1.8	1600	≥ 2.0	without spike, uncoated, with alligator teeth
JHY-FB-18-160-O-P-A1	1.8	1600	≥ 2.0	without spike, coated, with alligator teeth
JHY-FB-23-160-O-O-A1	2.3	1600	≥ 2.8	without spike, uncoated, with alligator teeth
JHY-FB-23-160-O-P-A1	2.3	1600	≥ 2.8	without spike, coated, with alligator teeth
JHY-FB-23-160-W-O-A1	2.3	1600	≥ 2.8	with spike, coated, with alligator teeth
JHY-FB-23-180-O-O-A1	2.3	1800	≥ 2.8	without spike, uncoated, with alligator teeth
JHY-FB-23-180-O-P-A1	2.3	1800	≥ 2.8	without spike, coated, with alligator teeth
JHY-FB-23-180-W-P-A1	2.3	1800	≥ 2.8	with spike, coated, with alligator teeth
JHY-FB-23-200-O-O-A1	2.3	2000	≥ 2.8	without spike, uncoated, with alligator teeth
JHY-FB-23-200-O-P-A1	2.3	2000	≥ 2.8	without spike, coated, with alligator teeth
JHY-FB-23-200-W-P-A1	2.3	2000	≥ 2.8	with spike, coated, with alligator teeth
JHY-FB-23-230-O-O-A1	2.3	2300	≥ 2.8	without spike, uncoated, with alligator teeth
JHY-FB-23-230-O-P-A1	2.3	2300	≥ 2.8	without spike, coated, with alligator teeth