

Disposable Grasping Forceps



Disposable Grasping Forceps



Forceps with alligator teeth

- Suit for grasping flat foreign body (e.g. coin)




Forceps with prongs

- Suit for grasping excised polyp, disease tissues and etc.



Forceps with rat teeth

- Suit for grasping foreign body in complex shape (e.g. retrieval of non-regular shaped stent)



Forceps with pelican teeth

- Suit for grasping special foreign body in big size (e.g. rubber, fabric)

Specifications

Model	O.D. (mm)	Length (mm)	Endoscope Channel (mm)	Characteristics
JHY-FG-18-70-A4	1.8	70	≥2.0	A4: Pelican Teeth
JHY-FG-23-160-A1	2.3	1600	≥2.8	A1: Alligator Teeth
JHY-FG-23-180-A1	2.3	1800	≥2.8	
JHY-FG-23-230-A1	2.3	2300	≥2.8	
JHY-FG-23-160-A2	2.3	1600	≥2.8	A2: Rat Teeth
JHY-FG-23-180-A2	2.3	1800	≥2.8	
JHY-FG-23-230-A2	2.3	2300	≥2.8	
JHY-FG-23-160-A3	2.3	1600	≥2.8	A3: Rat Teeth with Alligator Teeth Jaws
JHY-FG-23-180-A3	2.3	1800	≥2.8	
JHY-FG-23-230-A3	2.3	2300	≥2.8	
JHY-FG-23-160-A4	2.3	1600	≥2.8	A4: Pelican Teeth
JHY-FG-23-180-A4	2.3	1800	≥2.8	
JHY-FG-23-230-A4	2.3	2300	≥2.8	
JHY-FG-23-160-B1	2.3	1600	≥2.8	B1: 3 Round Prongs
JHY-FG-23-230-B1	2.3	2300	≥2.8	
JHY-FG-23-160-B2	2.3	1600	≥2.8	B2: 4 Round Prongs
JHY-FG-23-230-B2	2.3	2300	≥2.8	
JHY-FG-23-160-B3	2.3	1600	≥2.8	