

Disposable Grasping Forceps (With Net)



Disposable Grasping Forceps (With Net)



Intended use

- Grasp and remove foreign bodies from digestive tract.

Specifications

Model	O.D.(mm)	Length(mm)	Endoscope Channel(mm)	NET Size (mm)
JHY-FG-18-120-C1	1.8	1200	≥2.0	20*50
JHY-FG-18-120-C2	1.8	1200	≥2.0	25*55
JHY-FG-18-120-C3	1.8	1200	≥2.0	30*60
JHY-FG-18-160-C1	1.8	1600	≥2.0	20*50
JHY-FG-18-160-C2	1.8	1600	≥2.0	25*55
JHY-FG-18-160-C3	1.8	1600	≥2.0	30*60
JHY-FG-25-230-C1	2.5	2300	≥2.8	20*50
JHY-FG-25-230-C2	2.5	2300	≥2.8	25*55
JHY-FG-25-230-C3	2.5	2300	≥2.8	30*60
JHY-FG-25-230-C4	2.5	2300	≥2.8	40*75
JHY-FG-25-230-C5	2.5	2300	≥2.8	40*65