


Disposable Polypectomy Snare



Intended use

- Removal of polyps or other excess tissues in the digestive tract used with diathermic energy.

Features

- Made by high-strength medical braided wire, offering a precise & quick cutting properties, as well as an outstanding electric coagulation.
- Loop rotates synchronously by rotating the 3-ring handle, highly increase efficiency.
- Standard power connection, compatible with all main high-frequency devices on the market.
- Ergonomic design of 3-ring handle, easier to hold and use.

Specifications

Oval Snare

Model	O.D.(mm)	Length(mm)	Loop Width(mm)	Characteristics
JHY-SD-18-160-20-A1	1.8	1600	20	Standard Rotation
JHY-SD-23-160-10-A1	2.3	1600	10	
JHY-SD-23-230-10-A1	2.3	2300	10	
JHY-SD-23-160-15-A1	2.3	1600	15	
JHY-SD-23-230-15-A1	2.3	2300	15	
JHY-SD-23-160-20-A1	2.3	1600	20	
JHY-SD-23-230-20-A1	2.3	2300	20	
JHY-SD-23-160-25-A1	2.3	1600	25	
JHY-SD-23-230-25-A1	2.3	2300	25	
JHY-SD-23-160-30-A1	2.3	1600	30	
JHY-SD-23-230-30-A1	2.3	2300	30	
JHY-SD-23-160-35-A1	2.3	1600	35	
JHY-SD-23-230-35-A1	2.3	2300	35	
JHY-SD-23-160-50-A1	2.3	1600	50	
JHY-SD-23-230-50-A1	2.3	2300	50	
JHY-SD-23-230-20-A1R	2.3	2300	20	Graduated Rotation
JHY-SD-23-230-30-A1R	2.3	2300	30	

Crescent Snare

Model	O.D.(mm)	Length(mm)	Loop Width(mm)	Characteristics
JHY-SD-18-160-20-B1	1.8	1600	20	Standard Rotation
JHY-SD-23-160-10-B1	2.3	1600	10	
JHY-SD-23-230-10-B1	2.3	2300	10	
JHY-SD-23-160-15-B1	2.3	1600	15	
JHY-SD-23-230-15-B1	2.3	2300	15	
JHY-SD-23-160-20-B1	2.3	1600	20	
JHY-SD-23-230-20-B1	2.3	2300	20	
JHY-SD-23-160-25-B1	2.3	1600	25	
JHY-SD-23-230-25-B1	2.3	2300	25	
JHY-SD-23-160-30-B1	2.3	1600	30	
JHY-SD-23-230-30-B1	2.3	2300	30	

Hexagonal Snare

Model	O.D.(mm)	Length(mm)	Loop Width(mm)	Characteristics	
JHY-SD-18-160-20-C1	1.8	1600	20	Standard Rotation	
JHY-SD-23-160-10-C1	2.3	1600	10		
JHY-SD-23-230-10-C1	2.3	2300	10		
JHY-SD-23-160-15-C1	2.3	1600	15		
JHY-SD-23-230-15-C1	2.3	2300	15		
JHY-SD-23-160-20-C1	2.3	1600	20		
JHY-SD-23-230-20-C1	2.3	2300	20		
JHY-SD-23-160-25-C1	2.3	1600	25		
JHY-SD-23-230-25-C1	2.3	2300	25		
JHY-SD-23-160-30-C1	2.3	1600	30		
JHY-SD-23-230-30-C1	2.3	2300	30		
JHY-SD-23-230-20-C1R	2.3	2300	20		Graduated Rotation
JHY-SD-23-230-30-C1R	2.3	2300	30		