

Cleaning Brush


double ends cleaning brush

vapor cleaning brush

- Used for cleaning vapor channel

Intended use

- Cleaning of endoscope working and vapor channels



Specifications

Model	Diameter of Brush (mm)	Length (mm)	Endoscope Channel (mm)	Characteristics
JHY-CB-02-120-C	2	1200	2.0	Disposable, double ends, plastic shaft
JHY-CB-02-230-C	2	2300	2.0	
JHY-CB-03-120-C	3	1200	2.0	
JHY-CB-03-160-C	3	1600	2.0	
JHY-CB-03-230-C	3	2300	2.0	
JHY-CB-06-180-C	6	1800	2.8/3.2	
JHY-CB-06-230-C	6	2300	2.8/3.2	Disposable, single end, plastic shaft
JHY-CB-02-120-C1	2	1200	2.0	
JHY-CB-06-230-C1	6	2300	2.8/3.2	Reusable, Nylon brush, Plastic handle
JHY-CB-E	11/05	160	Vapor Channel	