



EF Rhino-Laryngoscope



- ◆ Handheld, Portable, Easy use
- ◆ Wireless Video Transmission
- ◆ No Optical Fiber



A **Authorised Representative**
IBS MEDICAL EQUIPMENTS SDN BHD
No. 7-G, 3rd Mile Square Business Centre,
No. 151, Old Klang Road, 58100 Wilayah
Persekutuan Kuala Lumpur.

T Tel : +603 7971 4221
W Whatsapp : +6012 346 3400

W Website : www.ibsmedic.com
E Email : sales@ibsmedic.com



EF Rhino-Laryngoscope



◆ **W**ireless video signal transmission

◆ **I**ntegrated unify display system. Portable for easy use



◆ **B**ig Working Channel

specifications	EF28	EF30	EF38	EF45
Field of View	120°	120°	120°	120°
Depth of View Field	3mm-200mm	3mm-200mm	3mm-200mm	3mm-200mm
Range of Tip Bending	(U) 180° (D) 130°	(U) 180° (D) 130°	(U) 180° (D) 130°	(U) 180° (D) 130°
Outer Diameter of Insertion Tube	Φ2.8mm	Φ3.0mm	Φ3.8mm	Φ4.5mm
Inner Diameter of Instrument Channel	/	Φ1.2mm	/	Φ2.2mm
Working Length	365mm	365mm	365mm	365mm

A **Authorised Representative**
IBS MEDICAL EQUIPMENTS SDN BHD
No. 7-G, 3rd Mile Square Business Centre,
No. 151, Old Klang Road, 58100 Wilayah
Persekutuan Kuala Lumpur.

T Tel : +603 7971 4221
W Whatsapp : +6012 346 3400

W Website : www.ibsmedic.com
E Email : sales@ibsmedic.com