





HD Flexible Digital BRONCHOSCOPE



Super Big Working Channel Φ 3.2mm



Applicable for bed-side treatment



Large Field of View



Specifications	VS100
Imaging resolution	1280x800
Video output	HDMI
Data storage	USB Disk, SD card(Built-in 8GB)
Inner/ Video format	JPG/ BMP/ MP4
Power	DC3.7V Li-ion Battery; 100~240V; 50/60Hz

Specifications	BF52H	BF58H	BF62H
Outer diameter of insertion tube	Φ 5.2mm	Φ 5.8mm	Φ 6.2mm
Inner diameter of instrument channel	Φ 2.2mm	Φ 2.9mm	Φ 3.2mm
Field of view	120°		
Depth of view field	3mm-200mm		
Range of tip Bending	(U)180° (D)130°		
Working length	600mm		