

Disposable Urological Guidewire

MDA Reg No. : GB2581824-170383

Ordering Information

Zebra Guide Wire

- ◆ Kink resistant Nickel – Titanium alloy core with precise advancement and reflection feeling
- ◆ Zebra stripes jacket for clearer movement observation
- ◆ Easily lubricated PTFE coating
- ◆ Pass the narrows easily

Ordering Information

Model	Description	Tip	O.D (in)	L (cm)
201CH2812	Ultra Smooth	5CM soft	0.028	120
201CH2815				150
201CH2818				180
201CH3212			0.032	120
201BM3215				150
201CH3218				180
201CH3512			0.035	120
201CH3515				150
201CH3518				180

Model	Description	Tip	O.D (in)	L (cm)
201BM2812	Zebra Wire	5CM soft	0.028	120
201BM2815				150
201BM2818				180
201BM3212			0.032	120
201BM3215				150
201BM3218				180
201BM3512			0.035	120
201BM3515				150
201BM3518				180

Ultra Smooth Guide Wire

- ◆ Hydrophilic and ultra soft distal tip design
- ◆ Enhanced visibility under fluoroscopy
- ◆ Hydrophilic and lubricious PVP coating helps the navigation



Dual Tips Guide Wire

- ◆ Dual soft tips
- ◆ Diverse application
- ◆ Enhance success

Hybrid Guide Wire

- ◆ Special processed tip, clearer under fluoroscopy and B-Scan
- ◆ Kink resistant nickel-titanium alloy core with precise advancement and reflection
- ◆ Zebra stripes could be observed clearly under endoscopes for movement tracking
- ◆ Soft and smooth hydrophilic tip passes the curve tract easily

Disposable Ureteral Stent

MDA Reg No. : GC3971325-195113

Ordering Information

High Efficiency Triangular Stent

- ◆ Patented stent, drain and prevent the new stones forming
- ◆ Clearly observed in X-rays
- ◆ Triangle Prism structure design and concave for better drainage space
- ◆ High retention pigtail shape tails reduce potential stent movement


Model	Type	O.D. (Fr)	Effective Length	Package
203B6DG	Triangle	6	25/27/29/31	No Guidewire
203B7DG	Triangle	7	25/27/29/31	No Guidewire
203A4DG	Classic	4	25/27/29/31	No Guidewire
203A5DG	Classic	4.8	25/27/29/31	No Guidewire
203A6DG	Classic	6	25/27/29/31	No Guidewire
203A7DG	Classic	7	25/27/29/31	No Guidewire
203A4TC	Classic	4	25/27/29/31	Coated Guidewire
203A4CH	Classic	4	25/27/29/31	Smooth Guidewire
203A4BM	Classic	4	25/27/29/31	Zebra Guidewire

Ordering Information

Model	Type	O.D. (Fr)	Effective Length	Package
203A5TC	Classic	4.8	25/27/29/31	Coated Guidewire
203A5CH	Classic	4.8	25/27/29/31	Smooth Guidewire
203A5BM	Classic	4.8	25/27/29/31	Zebra Guidewire
203A6TC	Classic	6	25/27/29/31	Coated Guidewire
203A6CH	Classic	6	25/27/29/31	Smooth Guidewire
203A6BM	Classic	6	25/27/29/31	Zebra Guidewire
203A7TC	Classic	7	25/27/29/31	Coated Guidewire
203A7CH	Classic	7	25/27/29/31	Smooth Guidewire
203A7BM	Classic	7	25/27/29/31	Zebra Guidewire

High Efficiency Triangular Stent

- ◆ Biocompatible polyurethane (PU) that softens > 35% at body temperature
- ◆ Larger inner lumen and thinner stent wall design to promote the drainage
- ◆ Cylindrical body matching to the ureter reduces the foreign body sensation




Disposable Ureteral Access Sheath

MDA Reg No. : GB2744723-157011

Negative Pressure Suction Sheath Patented Y Shape Connector

- ◆ Silicone sealing ring, prevent the liquid leakage
- ◆ Printed calibration for depth control
- ◆ Air ventilation Slot to regulate the pressure
- ◆ Bypass channel for stones and liquid drainage
- ◆ Discharges the stones by continuous infusion and negative pressure suction during the lithotripsy, stabilises the temperature to reduce irritation and injury to the tract and mucous



The flexible design of the dilator head allows for smooth passage through the ureteral stenosis area. The entire body of the dilatory and sheath tube can be clearly observed under X-ray, resulting in more accurate positioning



Flexional Tip Type

- ◆ >180° flexional with the ureteroscope
- ◆ Access renal pelvis and calices together with the flexible ureteroscope and regulate the kidney inner pressure
- ◆ Clean while absorbing, improving efficiency

Classic Type


- ◆ Smooth inside, not interference instrument
- ◆ Smooth coating, access tract easily
- ◆ Stiff and straight

Disposable PCNL Dilatation Kit

MDA Reg No. : GB4366724-174338


Regulatory Nephrostomy Kit

- ◆ Soft distal tip, less injury to the renal parenchyma, stiff guidewire, never lose the tract
- ◆ Calibration printed firmly on dilators' body for access depth observation. Visible under fluoroscopy and B-Scan.
- ◆ Specially designed dilators provide good handling feel, less friction force, and minimise tissue damage
- ◆ Lacerable sheath, could be teared and removed easily by pulling the two ears



A complete kit for diverse solutions

Negative Pressure Suction Type Nephrostomy Kit (Y Shape Connector)



Patented negative pressure suction sheath
Stone discharging and breaking simultaneously

- ◆ Air ventilation slot to regulate the pressure
- ◆ Silicon sealing ring, prevent the liquid leakage
- ◆ Bypass channel for stones and liquid drainage

Model	Description	I.D. (Fr)	L (cm)	Package
206JT14TZ	Negative pressure suction kit	14	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12 Fr dilators
206JT16TZ		16	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14 Fr dilators
206JT18TZ		18	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16 Fr dilators
206JT20TZ		20	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16/18 Fr dilators
206JT22TZ		22	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16/18/20 Fr dilators
206JT24TZ		24	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16/18/20/22 Fr dilators
206PT14TZ	Regular Kit	14	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12 Fr dilators
206PT16TZ		16	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14 Fr dilators
206PT18TZ		18	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16 Fr dilators
206PT20TZ		20	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16/18 Fr dilators
206PT22TZ		22	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16/18/20 Fr dilators
206PT24TZ		24	18	1 x puncture needle, 1 x guide wire, 1 x guide sheath, 1 x Y shape connector, 1 x 8/10/12/14/16/18/20/22 Fr dilators



Authorised Representative
IBS MEDICAL EQUIPMENTS SDN BHD
No. 7-G, 3rd Mile Square Business Centre,
No. 151, Old Klang Road, 58100 Wilayah
Persekutuan Kuala Lumpur.



Tel : +603 7971 4221



Whatsapp : +6012 346 3400



Website : www.ibsmedic.com



Email : sales@ibsmedic.com

Disposable Nitinol Stone Retrieval Basket

MDA Reg No. : GB10605424-175928

Features

◆ Memory Function


Shape memory Nickel-titanium alloy wire and high elasticity for easier stone retrieval

◆ Slim Design

The maximum width of the insertion part is 0.8mm, it passes the flexible urethroscope easily

◆ Reliability

High-performance materials and workmanship ensures the mesh basket's repeated opening and closing reliable without deformation



Two Options

- ◆ Three-Wire Basket: flat and opening head, easy way to capture and release large stones
- ◆ Four-Wire Basket: suitable to catch small stones from its lateral

Ordering Information

Model	Description	Basket Opening Width	Basket Length	Working Length	Maximum Width
211T10	Three-Wire	10mm	10mm	120cm	0.8
211F08	Four-Wire	8mm	18mm	120cm	0.8