

# RENT & OWN

- ✓ We offer rent and own of endoscopes
- ✓ Comprehensive coverage including replacement in 48 Hours
- ✓ Contract and coverage upto 48 Months



Cystoscope

Intubation scope / Bronchoscope

Rhionlaryngoscope

Ureteroscope

Video Endoscope	Model	Working length	Angulation	Photographing & Recording	Resolution	Contract (Months)	Outer diameter	Instrument channel
Flexible Video Bronchoscope	B41	600mm	U180°/D130°	Yes	400*400	24	3.2mm	1.2mm
							4.0mm	1.8mm
							5.2mm	2.6mm
Flexible Video Ureteroscope	U41	685mm	U275°/D275°	Yes	400*400	24	3.2mm	1.2mm
Flexible Video Cystoscope	C41	380mm	U210°/D130°	Yes	400*400	24	5.2mm	2.6mm
Video Portable Nasopharyngoscope	E30	330mm	U180°/D130°	Yes	400*400	24	2.8mm	/
							3.8mm	
Item	Model	Screen size	Storage	Photographing & Recording	Display Resolution	Contract (Months)	Battery	Continuous Operating Time
4.3" Monitor	M43	4.3" touch screen	8 G	Yes	800 X 480 (RGB)	24	3.7V lithium battery, rechargeable	>4 hours
10" Monitor	C100	10" touch screen	40 G	Yes	1280*800	24	Built-in Battery	>3 hours
Trolley	T1	5 Feet				24		

## Terms & Conditions

- ✓ Contract is 24 Months, no early termination of contract.
- ✓ All service will be covered by IBS Medic except for Mishandling cases
- ✓ In case of any issue, a loaner set will be deploy within 48 hours, the original equipments will be returned upon service
- ✓ Equipments will be owned by customer upon ending of contract



MDA Registration  
In progress

HQ Kuala Lumpur  
Unit 07-05, Level 07, Tower , Vertical Business  
Suite Avenue 3, Bangsar South No.8, Jalan  
Kerinci, 59200 Kuala Lumpur, Malaysia.  
T: 603 2242 4024

Penang Branch  
1-04-33 E-Gate, Lebu Tunku Kudin 2, 11700  
Gelugor, Penang, Malaysia.  
T: +604 824 9292