


Disposable Stone Extraction Basket



Disposable Stone Extraction Basket

Intended use

- Removal of gallstone in biliary duct and foreign bodies in upper and lower digestive tract

Features

- Innovative handle design, with functions of push, pull and rotation, easier to grasp gallstone and foreign body.
- Convenient for injection of contrast medium with injection port on handle.
- Made by advanced alloyed material, ensures good shape retention even after difficult stone removal.

Specifications

Model	O.D. (mm)	Length (mm)	Basket Width (mm)	Endoscope Channel (mm)	Characteristics
JHY-BAS-18-70-15-S4-A-O	1.8	700	15	>2.0	Stainless Steel, 4 wires, diamond shape, single-hand operation
JHY-BAS-18-70-20-S4-A-O	1.8	700	20	>2.0	
JHY-BAS-23-200-25-S4-A-D	2.3	2000	25	>2.8	Stainless Steel, 4 wires, diamond shape, both hands operation
JHY-BAS-23-200-30-S4-A-D	2.3	2000	30	>2.8	
JHY-BAS-23-200-35-S4-A-D	2.3	2000	35	>2.8	
JHY-BAS-18-70-15-N4-B-O	1.8	700	15	>2.0	Nitinol, 4 wires, spiral shape, single-hand operation
JHY-BAS-18-70-20-N4-B-O	1.8	700	20	>2.0	
JHY-BAS-23-200-20-N4-B-D	2.3	2000	20	>2.8	Nitinol, 4 wires, spiral shape, both hands operation
JHY-BAS-23-200-25-N4-B-D	2.3	2000	25	>2.8	
JHY-BAS-23-200-30-N4-B-D	2.3	2000	30	>2.8	
JHY-BAS-18-70-15-N4-C-O	1.8	700	15	>2.0	Nitinol, 4 wires, semi-circle shape, single-hand operation
JHY-BAS-18-70-20-N4-C-O	1.8	700	20	>2.0	
JHY-BAS-23-200-20-N4-C-D	2.3	2000	20	>2.8	Nitinol, 4 wires, semi-circle shape, both hands operation
JHY-BAS-23-200-25-N4-C-D	2.3	2000	25	>2.8	
JHY-BAS-23-200-30-N4-C-D	2.3	2000	30	>2.8	

