


Disposable Stone Extraction Balloon



Disposable Stone Extraction Balloon




Intended use

- Stone removal in biliary tract, including biliary sludge and residual gallstone left in biliary tract after Mechanical Lithotripsy.

Features

- Innovational design of three lumens (lumens for guide wire, injection and inflation) makes the insertion and radiography easier and more convenient.
- Three different sizes for each balloon, meets demands of different patients.
- Equipped with Radiopaque marks on both ends of balloon, convenient for locating under X-ray.
- Unique design of tapering catheter, improves the insertion into papilla.



Specifications

Model	O.D.(mm)	Length(mm)	Endoscope Channel(mm)	Balloon Size(mm)
JHY-BAL-23-8.5/11.5/15-B	2.3	2000	> 3.2	8.5/11.5/15
JHY-BAL-23-9/12/15-B	2.3	2000	> 3.2	9/12/15
JHY-BAL-23-11/13/15-B	2.3	2000	> 3.2	11/13/15
JHY-BAL-23-10/13/16-B	2.3	2000	> 3.2	10/13/16