

# Disposable Hot Biopsy Forceps



## Disposable Hot Biopsy Forceps

### Intended use

- Compatible with high-frequency surgical equipment and endoscope, it is used for peeling of tiny polyps or redundant tissues in digestive tract as well as for blood coagulation.



### Features

- Easy and comfortable to use.
- With smooth tube, easy to pass through the endoscope channel.
- With puncture-proof insulation, high safety.



### Specifications

Model	O.D. (mm)	Length (mm)	Endoscope Channel (mm)	Characteristics
JHY-FD-23-160-O	2.3	1600	≥2.8	Without spike
JHY-FD-23-180-O	2.3	1800	≥2.8	
JHY-FD-23-230-O	2.3	2300	≥2.8	
JHY-FD-23-260-O	2.3	2600	≥2.8	With spike
JHY-FD-23-160-W	2.3	1600	≥2.8	
JHY-FD-23-180-W	2.3	1800	≥2.8	
JHY-FD-23-230-W	2.3	2300	≥2.8	
JHY-FD-23-260-W	2.3	2600	≥2.8	